


24-Methylenecycloartane-3 β ,26-diol

Chemical Properties

CAS No. : 17020-27-8
 Formula: C₃₁H₅₂O₂
 Molecular Weight: 456.8
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	24-Methylenecycloartane-3 β ,26-diol is a natural product and can be used as a standard reference.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1891 mL	10.9457 mL	21.8914 mL
5 mM	0.4378 mL	2.1891 mL	4.3783 mL
10 mM	0.2189 mL	1.0946 mL	2.1891 mL
50 mM	0.0438 mL	0.2189 mL	0.4378 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481