


(±)-ε-Viniferin**Chemical Properties**

CAS No. :	253435-07-3
Formula:	C ₂₈ H ₂₂ O ₆
Molecular Weight:	454.47
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

**Biological Description**

Description	"(±)-ε-Viniferin, a racemate of ε-Viniferin, exhibits P450 inhibitory antioxidants, along with hepato-protective and antimicrobial activities. It also demonstrates antibacterial and antibiofilm activities against Gram-positive bacteria Streptococcus pneumoniae, with a minimum inhibitory concentration (MIC) of 20 μM [1]."
Targets(IC50)	Others,Antibacterial,Cytochromes P450

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2004 mL	11.0018 mL	22.0037 mL
5 mM	0.4401 mL	2.2004 mL	4.4007 mL
10 mM	0.220 mL	1.1002 mL	2.2004 mL
50 mM	0.044 mL	0.220 mL	0.4401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481