


Secotubeimoside I

Chemical Properties

CAS No. :	106235-32-9
Formula:	C ₆₀ H ₁₀₀ O ₃₀
Molecular Weight:	1337.5
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Secotubeimoside I is a tubeimoside I analogue and can be used as a standard reference.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7477 mL	3.7383 mL	7.4766 mL
5 mM	0.1495 mL	0.7477 mL	1.4953 mL
10 mM	0.0748 mL	0.3738 mL	0.7477 mL
50 mM	0.015 mL	0.0748 mL	0.1495 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481