

Isochlorogenic acid C 4'-O-glucoside

Chemical Properties

CAS No. :

Formula: C₃₁H₃₄O₁₇

Molecular Weight: 678.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Isochlorogenic acid C 4'-O-glucoside is a natural product that can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4736 mL	7.3681 mL	14.7362 mL
5 mM	0.2947 mL	1.4736 mL	2.9472 mL
10 mM	0.1474 mL	0.7368 mL	1.4736 mL
50 mM	0.0295 mL	0.1474 mL	0.2947 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481