

9,11-Dehydro-beta-boswellic acid

Chemical Properties


CAS No. : 471-65-8

Formula: C₃₀H₄₆O₃

Molecular Weight: 454.7

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	9,11-Dehydro-β-boswellic acid is a compound with urease inhibitory activity and acts as an inhibitor of 5-lipoxygenase (5-LO), a critical enzyme in leukotriene biosynthesis.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1993 mL	10.9963 mL	21.9925 mL
5 mM	0.4399 mL	2.1993 mL	4.3985 mL
10 mM	0.2199 mL	1.0996 mL	2.1993 mL
50 mM	0.044 mL	0.2199 mL	0.4399 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481