

(R)-5,3'-Dimethyl hesperidin

Chemical Properties

CAS No. :

Formula: C₃₀H₃₈O₁₅

Molecular Weight: 638.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Hesperidin analog serves as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5659 mL	7.8296 mL	15.6593 mL
5 mM	0.3132 mL	1.5659 mL	3.1319 mL
10 mM	0.1566 mL	0.783 mL	1.5659 mL
50 mM	0.0313 mL	0.1566 mL	0.3132 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481