


Hederagenin 3-O-(2-O-acetyl- α -L-arabinopyranoside)

Chemical Properties

CAS No. : 87562-05-8
 Formula: C₃₇H₅₈O₉
 Molecular Weight: 646.9
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hederagenin 3-O-(2-O-acetyl- α -L-arabinopyranoside) significantly elevated alkaline phosphatase (ALP) activity and enhanced the protein expression levels of bone sialoprotein (BSP) and osteocalcin (OC).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5458 mL	7.7292 mL	15.4583 mL
5 mM	0.3092 mL	1.5458 mL	3.0917 mL
10 mM	0.1546 mL	0.7729 mL	1.5458 mL
50 mM	0.0309 mL	0.1546 mL	0.3092 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481