


1-Phenylethyl beta-D-glucoside

Chemical Properties

CAS No. :	93199-03-2
Formula:	C ₁₄ H ₂₀ O ₆
Molecular Weight:	284.3
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Phenylethyl beta-D-glucoside is a natural product that can be used as a reference standard. The CAS number of 1-Phenylethyl beta-D-glucoside is 93199-03-2.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5174 mL	17.5871 mL	35.1741 mL
5 mM	0.7035 mL	3.5174 mL	7.0348 mL
10 mM	0.3517 mL	1.7587 mL	3.5174 mL
50 mM	0.0703 mL	0.3517 mL	0.7035 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481