


Ephedrannin D1

Chemical Properties

CAS No. :	1592431-55-4
Formula:	C ₃₀ H ₂₄ O ₁₃
Molecular Weight:	592.5
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ephedrannin D1 is a natural product that can be used as a reference standard. The CAS number of Ephedrannin D1 is 1592431-55-4.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6878 mL	8.4388 mL	16.8776 mL
5 mM	0.3376 mL	1.6878 mL	3.3755 mL
10 mM	0.1688 mL	0.8439 mL	1.6878 mL
50 mM	0.0338 mL	0.1688 mL	0.3376 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481