


Resinacein D

Chemical Properties

CAS No. :	2231061-47-3
Formula:	C ₃₀ H ₄₆ O ₅
Molecular Weight:	486.7
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Resinacein D is a natural product that can be used as a reference standard. The CAS number of Resinacein D is 2231061-47-3.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0547 mL	10.2733 mL	20.5465 mL
5 mM	0.4109 mL	2.0547 mL	4.1093 mL
10 mM	0.2055 mL	1.0273 mL	2.0547 mL
50 mM	0.0411 mL	0.2055 mL	0.4109 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481