

## Ganolucidic acid B

## Chemical Properties


CAS No. : 98683-75-1

Formula: C<sub>30</sub>H<sub>46</sub>O<sub>6</sub>

Molecular Weight: 502.7

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ganolucidic acids B demonstrated the ability to induce hPXR-mediated CYP3A4 expression.
Targets(IC50)	Others,Cytochromes P450

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9893 mL	9.9463 mL	19.8926 mL
5 mM	0.3979 mL	1.9893 mL	3.9785 mL
10 mM	0.1989 mL	0.9946 mL	1.9893 mL
50 mM	0.0398 mL	0.1989 mL	0.3979 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481