

3''-O-Caffeoylhesperidin

Chemical Properties

CAS No. :

Formula: C37H40O18

Molecular Weight: 772.7

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	3''-O-Caffeoylhesperidin is a natural product that can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2942 mL	6.4708 mL	12.9416 mL
5 mM	0.2588 mL	1.2942 mL	2.5883 mL
10 mM	0.1294 mL	0.6471 mL	1.2942 mL
50 mM	0.0259 mL	0.1294 mL	0.2588 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481