


Isorhamnetin 3,5-O-diglucoside

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 2035413-03-5 |
| Formula: | C ₂₈ H ₃₂ O ₁₇ |
| Molecular Weight: | 640.6 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|---------------|---|
| Description | Isorhamnetin 3,5-O-diglucoside is a natural product that can be used as a reference standard. The CAS number of Isorhamnetin 3,5-O-diglucoside is 2035413-03-5. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.561 mL | 7.8052 mL | 15.6104 mL |
| 5 mM | 0.3122 mL | 1.561 mL | 3.1221 mL |
| 10 mM | 0.1561 mL | 0.7805 mL | 1.561 mL |
| 50 mM | 0.0312 mL | 0.1561 mL | 0.3122 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481