


24-epi-7,8-Didehydrocimigenol 3-xyloside

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 150972-77-3 |
| Formula: | C ₃₅ H ₅₄ O ₉ |
| Molecular Weight: | 618.8 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|---------------|--|
| Description | 24-epi-7,8-Didehydrocimigenol 3-xyloside is a natural product that can be used as a reference standard. The CAS number of 24-epi-7,8-Didehydrocimigenol 3-xyloside is 150972-77-3. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.616 mL | 8.0802 mL | 16.1603 mL |
| 5 mM | 0.3232 mL | 1.616 mL | 3.2321 mL |
| 10 mM | 0.1616 mL | 0.808 mL | 1.616 mL |
| 50 mM | 0.0323 mL | 0.1616 mL | 0.3232 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481