


## Kihadanin D

## Chemical Properties

CAS No. :	2770024-83-2
Formula:	C <sub>27</sub> H <sub>32</sub> O <sub>9</sub>
Molecular Weight:	500.6
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Kihadanin D is a natural product that can be used as a reference standard. The CAS number of Kihadanin D is 2770024-83-2.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9976 mL	9.988 mL	19.976 mL
5 mM	0.3995 mL	1.9976 mL	3.9952 mL
10 mM	0.1998 mL	0.9988 mL	1.9976 mL
50 mM	0.040 mL	0.1998 mL	0.3995 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481