


Methyl ganoderate E

Chemical Properties

CAS No. :	98718-43-5
Formula:	C ₃₁ H ₄₂ O ₇
Molecular Weight:	526.7
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Methyl ganoderate E is a natural product that can be used as a reference standard. The CAS number of Methyl ganoderate E is 98718-43-5.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8986 mL	9.4931 mL	18.9861 mL
5 mM	0.3797 mL	1.8986 mL	3.7972 mL
10 mM	0.1899 mL	0.9493 mL	1.8986 mL
50 mM	0.038 mL	0.1899 mL	0.3797 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481