


Solafuranone

Chemical Properties

CAS No. :	367965-50-2
Formula:	C ₁₅ H ₂₀ O ₂
Molecular Weight:	232.3
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Solafuranone exhibits cytotoxic activity against human cancer cell lines at concentrations as high as 30 μ M and possesses anti-inflammatory properties.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3048 mL	21.5239 mL	43.0478 mL
5 mM	0.861 mL	4.3048 mL	8.6096 mL
10 mM	0.4305 mL	2.1524 mL	4.3048 mL
50 mM	0.0861 mL	0.4305 mL	0.861 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481