

Ganoderenic acid G

Chemical Properties


CAS No. : 120481-73-4

Formula: C₃₀H₄₀O₇

Molecular Weight: 512.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ganoderenic acid G is a natural product that can be used as a reference standard. The CAS number of Ganoderenic acid G is 120481-73-4.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9508 mL	9.7542 mL	19.5084 mL
5 mM	0.3902 mL	1.9508 mL	3.9017 mL
10 mM	0.1951 mL	0.9754 mL	1.9508 mL
50 mM	0.039 mL	0.1951 mL	0.3902 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481