


## Bayin

## Chemical Properties

CAS No. :	3681-96-7
Formula:	C <sub>21</sub> H <sub>20</sub> O <sub>9</sub>
Molecular Weight:	416.4
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Bayin is a natural product that can be used as a reference standard. The CAS number of Bayin is 3681-96-7.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4015 mL	12.0077 mL	24.0154 mL
5 mM	0.4803 mL	2.4015 mL	4.8031 mL
10 mM	0.2402 mL	1.2008 mL	2.4015 mL
50 mM	0.048 mL	0.2402 mL	0.4803 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481