


Schizotenuin A

Chemical Properties

CAS No. : 144608-09-3
 Formula: C₃₆H₂₈O₁₆
 Molecular Weight: 716.6
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Schizotenuin A is a natural product that can be used as a reference standard. The CAS number of Schizotenuin A is 144608-09-3.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3955 mL	6.9774 mL	13.9548 mL
5 mM	0.2791 mL	1.3955 mL	2.791 mL
10 mM	0.1395 mL	0.6977 mL	1.3955 mL
50 mM	0.0279 mL	0.1395 mL	0.2791 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481