


Ellagic acid 4-O-alpha-L-arabinofuranoside

Chemical Properties

CAS No. :	358617-39-7
Formula:	C ₁₉ H ₁₄ O ₁₂
Molecular Weight:	434.3
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ellagic acid 4-O-alpha-L-arabinofuranoside is a natural product that can be used as a reference standard. The CAS number of Ellagic acid 4-O-alpha-L-arabinofuranoside is 358617-39-7.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3026 mL	11.5128 mL	23.0256 mL
5 mM	0.4605 mL	2.3026 mL	4.6051 mL
10 mM	0.2303 mL	1.1513 mL	2.3026 mL
50 mM	0.0461 mL	0.2303 mL	0.4605 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481