


CDK-IN-11

## Chemical Properties

CAS No. : 2376876-87-6  
 Formula: C<sub>25</sub>H<sub>21</sub>BrN<sub>4</sub>O<sub>2</sub>  
 Molecular Weight: 489.36  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	CDK-IN-11, a heterocyclic compound, promotes cardiomyocyte maturation [1].
Targets(IC50)	Others,CDK

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0435 mL	10.2174 mL	20.4349 mL
5 mM	0.4087 mL	2.0435 mL	4.087 mL
10 mM	0.2043 mL	1.0217 mL	2.0435 mL
50 mM	0.0409 mL	0.2043 mL	0.4087 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481