

DPPE-MPEG(2000)

Chemical Properties

CAS No. : 459429-07-3

Formula: (C₂H₄O)_nC₃₉H₇₆N₁₀O₁₀P

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DPPE-MPEG(2000) is a polyethylene glycol (PEG) derivative of 1,2-dipalmitoyl-rac-glycero-3-phosphoethanolamine (DPPE), characterized as a PEGylated phospholipid.
Targets(IC50)	Others,Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481