

## Hydroxy Lenalidomide

## Chemical Properties


CAS No. : 1421593-78-3

Formula: C<sub>13</sub>H<sub>13</sub>N<sub>3</sub>O<sub>4</sub>

Molecular Weight: 275.26

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Hydroxy lenalidomide, a metabolite of the thalidomide analog and anticancer agent lenalidomide, is a derivative substance involved in the metabolic process of its parent compound.
Targets(IC50)	Others, Drug Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6329 mL	18.1646 mL	36.3293 mL
5 mM	0.7266 mL	3.6329 mL	7.2659 mL
10 mM	0.3633 mL	1.8165 mL	3.6329 mL
50 mM	0.0727 mL	0.3633 mL	0.7266 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481