

C26:1 Sphingomyelin (d18:1/26:1(17Z))

Chemical Properties


CAS No. : 1448012-56-3

Formula: C₄₉H₉₇N₂O₆P

Molecular Weight: 841.3

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	C26:1 Sphingomyelin (d18:1/26:1(17Z)) is a lipid molecule that can be used in life science related research. The CAS number of C26:1 Sphingomyelin (d18:1/26:1(17Z)) is 1448012-56-3.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1886 mL	5.9432 mL	11.8864 mL
5 mM	0.2377 mL	1.1886 mL	2.3773 mL
10 mM	0.1189 mL	0.5943 mL	1.1886 mL
50 mM	0.0238 mL	0.1189 mL	0.2377 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481