


## 1-Palmitoyl-2-Oleoyl-3-Linolenoyl-rac-glycerol

## Chemical Properties

CAS No. :	78144-20-4
Formula:	C55H98O6
Molecular Weight:	855.36
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Palmitoyl-2-oleoyl-3-linolenoyl-rac-glycerol is a triacylglycerol comprising palmitic acid, oleic acid, and $\alpha$ -linolenic acid esterified at the sn-1, sn-2, and sn-3 positions, respectively.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1691 mL	5.8455 mL	11.691 mL
5 mM	0.2338 mL	1.1691 mL	2.3382 mL
10 mM	0.1169 mL	0.5845 mL	1.1691 mL
50 mM	0.0234 mL	0.1169 mL	0.2338 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481