


2-Myristyldistearin

Chemical Properties

CAS No. : 57396-99-3
Formula: C₅₃H₁₀₂O₆
Molecular Weight: 835.37
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Myristyldistearin (SMS), a triacylglycerol comprising myristic acid and stearic acid [1], showcases its intricate structural composition.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1971 mL	5.9854 mL	11.9707 mL
5 mM	0.2394 mL	1.1971 mL	2.3941 mL
10 mM	0.1197 mL	0.5985 mL	1.1971 mL
50 mM	0.0239 mL	0.1197 mL	0.2394 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481