


## 1-Arachidoyl-rac-glycerol

## Chemical Properties

CAS No. :	50906-68-8
Formula:	C23H46O4
Molecular Weight:	386.61
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Arachidoyl-rac-glycerol, a monoacylglycerol, features arachidic acid at the sn-1 position and has been detected in pork muscle.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5866 mL	12.9329 mL	25.8659 mL
5 mM	0.5173 mL	2.5866 mL	5.1732 mL
10 mM	0.2587 mL	1.2933 mL	2.5866 mL
50 mM	0.0517 mL	0.2587 mL	0.5173 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481