

## 1,2-Diarachidoyl-rac-glycerol

## Chemical Properties


CAS No. : 89648-24-8

Formula: C43H84O5

Molecular Weight: 681.12

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1,2-Diarachidoyl-rac-glycerol is a diacylglycerol featuring arachidic acid at both the sn-1 and sn-2 positions.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4682 mL	7.3409 mL	14.6817 mL
5 mM	0.2936 mL	1.4682 mL	2.9363 mL
10 mM	0.1468 mL	0.7341 mL	1.4682 mL
50 mM	0.0294 mL	0.1468 mL	0.2936 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481