


## 1,2-Dielaidoyl-3-Stearoyl-rac-glycerol

## Chemical Properties

CAS No. :	93452-41-6
Formula:	C57H106O6
Molecular Weight:	887.5
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1,2-Dielaidoyl-3-stearoyl-rac-glycerol is a triacylglycerol comprised of elaidic acid at the (sn-1) and (sn-2) positions, and stearic acid at the (sn-3) position.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1268 mL	5.6338 mL	11.2676 mL
5 mM	0.2254 mL	1.1268 mL	2.2535 mL
10 mM	0.1127 mL	0.5634 mL	1.1268 mL
50 mM	0.0225 mL	0.1127 mL	0.2254 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481