


1,2-Dilinoleoyl-3- γ -Linolenoyl-rac-glycerol

Chemical Properties

CAS No. :	852369-68-7
Formula:	C57H96O6
Molecular Weight:	877.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,2-Dilinoleoyl-3- γ -linolenoyl-rac-glycerol is a triacylglycerol featuring linoleic acid at its sn-1 and sn-2 positions and γ -linolenic acid at the sn-3 position.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1398 mL	5.6988 mL	11.3977 mL
5 mM	0.228 mL	1.1398 mL	2.2795 mL
10 mM	0.114 mL	0.5699 mL	1.1398 mL
50 mM	0.0228 mL	0.114 mL	0.228 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481