

## DSPE-MPEG(2000) sodium

## Chemical Properties

CAS No. : 459428-35-4  
Formula: (C<sub>2</sub>H<sub>4</sub>O)<sub>n</sub>C<sub>43</sub>H<sub>83</sub>N<sub>0</sub>O<sub>10</sub>P?Na  
Molecular Weight:  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-MPEG(2000), a PEGylated derivative of 1,2-distearoyl-sn-glycero-3-PE (DSPE), is employed in the formulation of liposomes for delivering anticancer and antimalarial agents.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481