


Varinolic Acid

Chemical Properties

CAS No. :	4707-50-0
Formula:	C ₁₀ H ₁₂ O ₄
Molecular Weight:	196.1999
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Varinolic acid, an analytical reference standard, plays a crucial role as an intermediate in the phytocannabinoid biosynthetic pathway. It serves as a precursor for synthesizing cannabidivarin (CBDV), cannabigerovarin (CBGV), and cannabivarin (CBV). This compound is designed for research and forensic applications.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0968 mL	25.4842 mL	50.9684 mL
5 mM	1.0194 mL	5.0968 mL	10.1937 mL
10 mM	0.5097 mL	2.5484 mL	5.0968 mL
50 mM	0.1019 mL	0.5097 mL	1.0194 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481