


## KT109 N2 Regioisomer

## Chemical Properties

CAS No. :	2749638-90-0
Formula:	C <sub>27</sub> H <sub>26</sub> N <sub>4</sub> O
Molecular Weight:	422.52
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	KT109 N2 regioisomer, an N2-carbamoylated variant of the diacylglycerol lipase $\beta$ (DAGL $\beta$ ) inhibitor KT109, serves as an internal standard for KT109 quantification.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3668 mL	11.8338 mL	23.6675 mL
5 mM	0.4734 mL	2.3668 mL	4.7335 mL
10 mM	0.2367 mL	1.1834 mL	2.3668 mL
50 mM	0.0473 mL	0.2367 mL	0.4734 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481