


## Piperitenone

## Chemical Properties

CAS No. :	491-09-8
Formula:	C <sub>10</sub> H <sub>14</sub> O
Molecular Weight:	150.2176
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Piperitenone, a monoterpene identified in various plants, exhibits antioxidant activity and inhibits the peroxidation of linoleic acid.
Targets(IC50)	Others, Mitochondrial Metabolism, ROS

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.6569 mL	33.2845 mL	66.569 mL
5 mM	1.3314 mL	6.6569 mL	13.3138 mL
10 mM	0.6657 mL	3.3285 mL	6.6569 mL
50 mM	0.1331 mL	0.6657 mL	1.3314 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481