

Compound N039-0045

Chemical Properties


CAS No. :

Formula: C₄₁H₆₈O₁₄

Molecular Weight: 784.981

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N039-0045 is a useful organic compound for research related to life sciences and the catalog number is T131265.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2739 mL	6.3696 mL	12.7392 mL
5 mM	0.2548 mL	1.2739 mL	2.5478 mL
10 mM	0.1274 mL	0.637 mL	1.2739 mL
50 mM	0.0255 mL	0.1274 mL	0.2548 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481