

Haloperidol-D4 trifluoroacetate

Chemical Properties


CAS No. :

Formula: C₂₃H₂₀D₄ClF₄NO₄

Molecular Weight: 493.912

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Haloperidol-D4 trifluoroacetate is a deuterated compound of Haloperidol trifluoroacetic acid.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0247 mL	10.1233 mL	20.2466 mL
5 mM	0.4049 mL	2.0247 mL	4.0493 mL
10 mM	0.2025 mL	1.0123 mL	2.0247 mL
50 mM	0.0405 mL	0.2025 mL	0.4049 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481