

## Anti-IFNA17 Monoclonal Antibody

### Product Details

Ig Type:	Rabbit monoclonal IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	150 kDa
Purification:	Protein A Affinity Purified

### Applications

Verified Activity:	Flow cytometry analysis of IFNA17 overexpressed 376F cells with TMAZ-0025, followed by goat anti-rabbit IgG-ABflo 647 (red line). The isotype control is rabbit IgG (black line).
Application:	ELISA,FCM,IF
Recommended	0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA

### Properties

Purity:	> 95% as determined by SDS-PAGE.
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	IFNA17
Antigen Species:	Human
Gene ID:	3451
Uniprot ID:	P01571
Synonyms:	Interferon alpha-T;Interferon alpha-I';Interferon alpha-17;Interferon alpha-88
Biology Area:	Immunology Research

### Research Background

Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481