

## Anti-CD200 Monoclonal Antibody

### Product Details

Molecular Weight: 150 kDa

### Properties

### Antigen Details

Synonyms: OX-2 membrane glycoprotein

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481

## Anti-CD200 Monoclonal Antibody

### Product Details

Ig Type:	Rabbit monoclonal IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	150 kDa
Purification:	Protein A Affinity Purified

### Applications

Verified Activity:	Flow cytometry analysis of CD200 overexpressed 293F cells with TMAZ-0047BI, which was counjugated with SA-PE (red). The WT 293F control is black line.
Application:	ELISA,FCM
Recommended	0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA

### Properties

Purity:	> 95% as determined by SDS-PAGE.
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	CD200
Antigen Species:	Human
Gene ID:	4345
Uniprot ID:	P41217
Synonyms:	OX-2 membrane glycoprotein
Biology Area:	Immunology Research

### Research Background

Costimulates T-cell proliferation. May regulate myeloid cell activity in a variety of tissues

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481