

Anti-AOAH Monoclonal Antibody

Product Details

Molecular Weight: 150 kDa

Properties

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481

Anti-AOAH Monoclonal Antibody

Product Details

| | |
|-------------------|-----------------------|
| Ig Type: | Rabbit monoclonal IgG |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Molecular Weight: | 150 kDa |

Applications

| | |
|--------------|--|
| Application: | FCM |
| Recommended | 0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA |

Properties

| | |
|----------------------|---|
| Purity: | > 95% as determined by SDS-PAGE. |
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|---------------------|
| Immunogen: | AOAH |
| Antigen Species: | Human |
| Gene ID: | 313 |
| Biology Area: | Immunology Research |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481