

Anti-DR5 Monoclonal Antibody

Product Details

Ig Type:	Rabbit monoclonal IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	150 kDa
Purification:	Protein A Affinity Purified

Applications

Verified Activity:	Flow cytometry analysis of DR5 overexpressed 293F cells with TMAZ-0134BI, followed by goat anti-rabbit IgG-ABflo 647 (red). The isotype control is rabbit IgG (black line).
Application:	ELISA,FCM,IF
Recommended	0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA

Properties

Purity:	> 95% as determined by SDS-PAGE.
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	DR5
Antigen Species:	Human
Gene ID:	8795
Uniprot ID:	O14763
Synonyms:	TNF-related apoptosis-inducing ligand receptor 2;Tumor necrosis factor receptor superfamily member 10B;Death receptor 5
Biology Area:	Cancer Research

Research Background

Receptor for the cytotoxic ligand TNFSF10/TRAIL (PubMed: 10549288). The adapter molecule FADD recruits caspase-8 to the activated receptor. The resulting death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiates the subsequent cascade of caspases (aspartate-specific cysteine proteases) mediating apoptosis. Promotes the activation of NF-kappa-B. Essential for ER stress-induced apoptosis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481