

Anti-CD44 Monoclonal Antibody

Product Details

Ig Type:	Rabbit monoclonal IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	150 kDa
Purification:	Protein A Affinity Purified

Applications

Verified Activity:	Flow cytometry analysis of CD44 overexpressed 367F cells with TMAZ-0172, followed by goat anti-rabbit IgG-ABflo 647 (red line). The isotype control is rabbit IgG (black line).
Application:	ELISA, FCM
Recommended	0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA

Properties

Purity:	> 95% as determined by SDS-PAGE.
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	CD44
Antigen Species:	Human
Gene ID:	960
Uniprot ID:	P16070-4
Synonyms:	Hermes antigen;Phagocytic glycoprotein I;Isoform 4 of CD44 antigen;Phagocytic glycoprotein 1;GP90 lymphocyte homing/adhesion receptor;Hyaluronate receptor;CDw44;HUTCH-I;Heparan sulfate proteoglycan;Epican;Extracellular matrix receptor III
Biology Area:	Drug Target

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481

Anti-CD44 Monoclonal Antibody

Product Details

Molecular Weight: 150 kDa

Properties

Antigen Details

Synonyms:

Hermes antigen;Phagocytic glycoprotein I;Isoform 4 of CD44 antigen;Phagocytic glycoprotein 1;GP90 lymphocyte homing/adhesion receptor;Hyaluronate receptor;CDw44;HUTCH-I;Heparan sulfate proteoglycan;Epican;Extracellular matrix receptor III

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481