

Anti-TGOLN2 Monoclonal Antibody

Product Details

Ig Type:	Rabbit monoclonal IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	150 kDa
Purification:	Protein A Affinity Purified

Applications

Verified Activity:	Flow cytometry analysis of TGOLN2 overexpressed 347F cells with TMAZ-0183, followed by goat anti-rabbit IgG-ABflo 647 (red line). The isotype control is rabbit IgG (black line).
Application:	ELISA, FCM, IF
Recommended	0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA

Properties

Purity:	> 95% as determined by SDS-PAGE.
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	TGOLN2
Antigen Species:	Human
Gene ID:	10618
Uniprot ID:	O43493
Synonyms:	Trans-Golgi network protein 2;Trans-Golgi network glycoprotein 48;Trans-Golgi network integral membrane protein 2;Trans-Golgi network glycoprotein 51;Trans-Golgi network glycoprotein 46
Biology Area:	Drug Target

Research Background

May be involved in regulating membrane traffic to and from trans-Golgi network

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481