

## Anti-TFF1 Monoclonal Antibody

### Product Details

Ig Type:	Rabbit monoclonal IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	150 kDa
Purification:	Protein A Affinity Purified

### Applications

Verified Activity:	Flow cytometry analysis of TFF1 overexpressed 407F cells with TMAZ-0214, followed by goat anti-rabbit IgG-ABflo 647 (red line). The isotype control is rabbit IgG (black line).
Application:	ELISA,FCM
Recommended	0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA

### Properties

Purity:	> 95% as determined by SDS-PAGE.
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	TFF1
Antigen Species:	Human
Gene ID:	7031
Uniprot ID:	P04155
Synonyms:	Trefoil factor 1;Polypeptide P1.A;PNR-2;Breast cancer estrogen-inducible protein;Protein pS2
Biology Area:	Cytokines

### Research Background

Stabilizer of the mucous gel overlying the gastrointestinal mucosa that provides a physical barrier against various noxious agents. May inhibit the growth of calcium oxalate crystals in urine

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481