

Anti-KERA Monoclonal Antibody

## Product Details

Molecular Weight: 150 kDa

## Properties

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481

## Anti-KERA Monoclonal Antibody

### Product Details

|                   |                       |
|-------------------|-----------------------|
| Ig Type:          | Rabbit monoclonal IgG |
| Reactivity:       | Human                 |
| Conjugation:      | Unconjugated          |
| Molecular Weight: | 150 kDa               |

### Applications

|              |                                                  |
|--------------|--------------------------------------------------|
| Application: | FCM                                              |
| Recommended  | 0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA |

### Properties

|                      |                                                                           |
|----------------------|---------------------------------------------------------------------------|
| Purity:              | > 95% as determined by SDS-PAGE.                                          |
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping:            | Shipping with blue ice.                                                   |

### Antigen Details

|                  |                 |
|------------------|-----------------|
| Immunogen:       | KERA            |
| Antigen Species: | Human           |
| Gene ID:         | 11081           |
| Biology Area:    | Other Functions |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481