

Anti-ATRAID (ECD 169AA) Monoclonal Antibody-Biotin

Product Details

Molecular Weight: 150 kDa

Properties

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481

Anti-ATRAID (ECD 169AA) Monoclonal Antibody-Biotin

Product Details

Ig Type: Rabbit monoclonal IgG
Reactivity: Human
Conjugation: Biotin
Molecular Weight: 150 kDa

Applications

Application: FCM
Recommended: 0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA

Properties

Purity: > 95% as determined by SDS-PAGE.
Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: ATRAID (ECD 169AA)
Antigen Species: Human
Gene ID: 51374
Biology Area: Immunology Research

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481