

## Anti-TPBG Monoclonal Antibody-Biotin

### Product Details

Molecular Weight: 150 kDa

### Properties

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481

## Anti-TPBG Monoclonal Antibody-Biotin

### Product Details

|                   |                             |
|-------------------|-----------------------------|
| Ig Type:          | Rabbit monoclonal IgG       |
| Reactivity:       | Human                       |
| Conjugation:      | Biotin                      |
| Molecular Weight: | 150 kDa                     |
| Purification:     | Protein A Affinity Purified |

### Applications

|                    |  |
|--------------------|--|
| Verified Activity: | Flow cytometry analysis of TPBG overexpressed 374F cells with TMAZ-0115B, followed by goat anti-rabbit IgG-ABflo 647 (red line). The isotype control is rabbit IgG (black line). |
| Application:       | ELISA,FCM  |
| Recommended        | 0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA   |

### Properties

|                      |   |
|----------------------|---|
| Purity:              | > 95% as determined by SDS-PAGE.  |
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping:            | Shipping with blue ice.   |

### Antigen Details

|                  |                 |
|------------------|-----------------|
| Immunogen:       | TPBG            |
| Antigen Species: | Human           |
| Gene ID:         | 7162            |
| Biology Area:    | Cancer Research |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481