

## Anti-hIGF2 Monoclonal Antibody-Biotin

### Product Details

Ig Type: Rabbit monoclonal IgG  
Reactivity: Human  
Conjugation: Biotin  
Molecular Weight: 150 kDa

### Applications

Application: FCM  
Recommended: 0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA

### Properties

Purity: > 95% as determined by SDS-PAGE.  
Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: hIGF2  
Antigen Species: Human  
Gene ID: 3481  
Biology Area: Immunology Research

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481

Anti-hIGF2 Monoclonal Antibody-Biotin

## Product Details

Molecular Weight: 150 kDa

## Properties

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481