

## Anti-MUCL1 Monoclonal Antibody-Biotin

### Product Details

Molecular Weight: 150 kDa

### Properties

### Antigen Details

Synonyms: Mucin-like protein 1;Protein BS106;Small breast epithelial mucin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481

## Anti-MUCL1 Monoclonal Antibody-Biotin

### Product Details

Ig Type:	Rabbit monoclonal IgG
Reactivity:	Human
Conjugation:	Biotin
Molecular Weight:	150 kDa
Purification:	Protein A Affinity Purified

### Applications

Verified Activity:	Flow cytometry analysis of MUCL1 overexpressed 474F cells with TMAZ-0266B, followed by goat anti-rabbit IgG-ABflo 647 (red line). The isotype control is rabbit IgG (black line).
Application:	ELISA, FCM
Recommended	0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA

### Properties

Purity:	> 95% as determined by SDS-PAGE.
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	MUCL1
Antigen Species:	Human
Gene ID:	118430
Uniprot ID:	Q96DR8
Synonyms:	Mucin-like protein 1;Protein BS106;Small breast epithelial mucin
Biology Area:	Other Functions

### Research Background

May play a role as marker for the diagnosis of metastatic breast cancer

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481