

Anti-JAM2 Monoclonal Antibody-Biotin

Product Details

| | |
|-------------------|-----------------------------|
| Ig Type: | Rabbit monoclonal IgG |
| Reactivity: | Human |
| Conjugation: | Biotin |
| Molecular Weight: | 150 kDa |
| Purification: | Protein A Affinity Purified |

Applications

| | |
|--------------------|--|
| Verified Activity: | Flow cytometry analysis of JAM2 overexpressed 467F cells with TMAZ-0270B, followed by goat anti-rabbit IgG-ABflo 647 (red line). The isotype control is rabbit IgG (black line). |
| Application: | ELISA,FCM |
| Recommended | 0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA |

Properties

| | |
|----------------------|---|
| Purity: | > 95% as determined by SDS-PAGE. |
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|---|
| Immunogen: | JAM2 |
| Antigen Species: | Human |
| Gene ID: | 58494 |
| Uniprot ID: | P57087 |
| Synonyms: | Junctional adhesion molecule 2;Junctional adhesion molecule B;Vascular endothelial junction-associated molecule |
| Biology Area: | Other Functions |

Research Background

Junctional adhesion protein that mediates heterotypic cell-cell interactions with its cognate receptor JAM3 to regulate different cellular processes (PubMed: 11590146, PubMed: 11823489, PubMed: 24357068). Plays a role in homing and mobilization of hematopoietic stem and progenitor cells within the bone marrow (PubMed: 24357068). At the surface of bone marrow stromal cells, it contributes to the retention of the hematopoietic stem and progenitor cells expressing JAM3 (PubMed: 11590146, PubMed: 24357068). Plays a central role in leukocytes extravasation by facilitating not only transmigration but also tethering and rolling of leukocytes along the endothelium (PubMed: 12239159). Tethering and rolling of leukocytes are dependent on the binding by JAM2 of the integrin alpha-4/beta-1 (PubMed: 12070135). Plays a role in spermatogenesis where JAM2 and JAM3, which are respectively expressed by Sertoli and germ cells, mediate an interaction between both cell types and play an essential role in the anchorage of germ cells onto Sertoli cells and the assembly of cell polarity complexes during spermatid differentiation (By similarity). Also functions as an inhibitory somatodendritic cue that prevents the myelination of non-axonal parts of neurons (By similarity). During myogenesis, it is involved in myocyte fusion (By similarity). May also play a role in

angiogenesis (By similarity)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481